

Product Name:	FRF302WWW/0GP
Description:	<p>Corepoint Scientific™ General Purpose Flammable Storage Refrigerator/Freezer Combination 30 Cu. Ft. Capacity (25 Ref./5 Freezer) One solid swing door/one interior solid swing door-freezer chamber only, right hinged, not reversible Mechanical thermostat Adjustable operating temperature range: Refrigerator: 1°C – 10°C Freezer: -15°C – -25°C Three adjustable refrigerator shelves Keyed door lock Magnetic door gasket for positive seal Stainless steel door handle Leveling legs Probe access port (3/8") White textured cabinet and door High density urethane foam cabinet and door insulation CFC-free refrigerant (R404a), foam insulation and packaging One year parts and labor warranty Exterior dimensions: 35 ¾"W x 35 ½"D x 72 ½"H Manual defrost 115V, 60Hz, 6 Amps, 1/3 HP UL/C-UL listed Units meet NFPA guidelines 45 and 70 Shipping weight: 396 lbs.</p>