

Product Name:	GPR723SSS/0
Description:	<p>Corepoint Scientific™ General Purpose Stainless-Steel Refrigerator  72 Cu. Ft. capacity  Three swing solid doors, self-closing, not reversible  9 ea. Adjustable Shelves  Digital display microprocessor temperature controller  Adjustable operating temperature: 1°C to 10°C (34°F to 50°F)  Audible and visual Hi/Lo temperature alarm  2 probes; 1 in air and 1 in sample bottle  Sensor error alarm  Continuous product (simulated) temperature display  Remote alarm contacts  Forced draft air circulation  Keyed door locks  Magnetic door gaskets for positive seal  Integrated door handle  LED interior lights are safety shielded and switch controlled  Probe access port (3/4")  Stainless steel exterior  White-powder coated aluminum interior with 304 stainless steel floor  Easily cleaned smooth interior  Hydrocarbon natural refrigerant (R290)  2 year parts and labor warranty, plus an additional 3 year  compressor parts warranty  6 swivel casters (5 3/4")  Exterior dimensions: 81" W x 33 3/4" D x 83 1/4" H with casters; 77  1/2" H without casters  Cycle defrost  115V, 60 Hz, 7.5 Amps, 1/2 HP  ETL listed  Shipping weight: 666 lbs.</p>