

Product Name:	LF301WWW/0M
Description:	<p>Corepoint Scientific™ General Purpose Manual Defrost Laboratory Freezer 30 Cu. Ft. Capacity One solid swing door, right hinged, not reversible Mechanical thermostat Adjustable operating temperature range: -15°-25°C Four fixed shelves Keyed door lock Magnetic door gasket for positive seal Stainless steel door handle Leveling legs Probe access port (3/8") White textured cabinet and door High density urethane foam cabinet and door insulation CFC-free refrigerant (R404a), foam insulation and packaging One year parts and labor warranty; plus additional four year compressor parts warranty Exterior dimensions: 35 3/4" W x 34 1/4" D x 72 1/2" H Manual defrost 115V, 60Hz, 6 Amps, 1/3 HP UL/C-UL listed Shipping weight: 365 lbs.</p>