

Product Name:	LR301WWW/0
Description:	<p>Corepoint Scientific™ General Purpose Refrigerator 30 Cu. Ft. Capacity One solid swing door, right hinged, not reversible Mechanical thermostat Temperature range: 0°C to 10°C Five adjustable shelves Keyed door lock Magnetic door gasket for positive seal Stainless steel door handle Leveling legs Probe access port (3/8") Powder coat paint exterior High density urethane foam cabinet and door insulation CFC-free refrigerant (R134a), foam insulation and packaging Two year parts and labor warranty, plus an additional three year compressor parts warranty Exterior dimensions: 35 3/4"W x 35" D x 72 1/2"H Auto defrost 115V, 60Hz, 6.8 amps, 1/3 HP UL/C-UL listed Shipping weight: 396 lbs.</p>