

Product Name:	LRF201WWW/0
Description:	<p>Corepoint Scientific™ General Purpose Refrigerator & Freezer Combination 20 cu. ft. capacity (15.5 ref / 5.5 freezer) Two exterior solid swing doors, right hinged, field reversible Mechanical thermostat Refrigerator temperature range : 0°C to 10°C Freezer chamber: -15°C to -10°C with two auto defrost cycles warming to 0°C per day Three shelves, refrigerator (one fixed, two adjustable); one removable shelf, freezer Magnetic door gasket for positive seal Heavy duty door handles Leveling legs Interior refrigerator light Two probe access ports (3/8"); one, refrigerator/one, freezer White textured cabinet and door High density urethane foam cabinet and door insulation CFC-free refrigerant (R134a), foam insulation and packaging Two year parts and labor warranty, plus an additional three year compressor parts warranty Exterior dimensions: 30"W x 34 ¼"D x 68 ½ " H Auto defrost 115V, 60 Hz, 6 Amps, ¼ HP Shipping weight: 231 lbs.</p>