<b>F</b>	
Product Name:	NSRI231WSW/0H
Description:	Corepoint Scientific™ Humidity and Temperature Stability Chamber 23 Cu. Ft. Capacity
	One swing solid door, self-closing, right hinged, not reversible Three,
	304 stainless steel sliding tray shelves
	Controller Information:
	Proportional integral derivative (PID) controller with LCD display
	Programmable ramp and soak
	Process phase graphic animations
	USB data port
	RS-485 communication (MODBUS-RTU)
	Real time clock
	Battery backup
	Event logs: 100
	Visual alarm indicator
	Audible alarm indicator
	Process phase specific temp alarm
	Process phase specific humidity alarm
	Power failure alarm
	Sensor failure alarm
	Temperature setpoint range: 4°C to 70°C
	Temperature control: +/- 0.1°C
	Temperature stability: +/- 0.5°C from 4°C to 70°C (Range
	dependent per test protocol)
	Temperature uniformity: +/- 1°C from 4°C to 70°C (Range
	dependent per test protocol)
	Humidity set point range: 15 to 95% RH
	Humidity control: +/- 0.1%RH
	Humidity Stability: +/- 3% RH (Range dependent per test protocol)
	Cabinet is foamed-in-place with EPA compliant high density;
	polyurethane foam insulation White-powder coated steel exterior
	304 stainless steel interior
	Plenum air distribution
	Floor drain
	Top mounted, variable speed compressor
	Hydrocarbon natural refrigerant
	Ultrasonic humidification
	Multi-element, independent proportional electric heating
	Temp sensor Type: NTC
	RH sensor Type: Capacitive
	Condensate evaporation
	Shielded, switched interior LED lighting
	Heavy duty door handle
	Magnetic door gasket for positive seal

Keyed door lock

Probe access port (3/4")

4 swivel casters (3 1/2")

18 months parts and labor warranty

3 year compressor warranty

Exterior dimensions: 26 7/8" W x 36 3/8" D x 80 3/4" H with casters;

77 1/4" H without casters

Interior dimensions: 22" W x 25 1/2" D x 59" H 10-60 PSI demineralized water supply required

115V, 60 Hz, 12.4 Amps, 1/3 HP

ETL Listed

Shipping weight: 380 lbs.