

<b>Product Name:</b>	NSWDR492WWG/0
<b>Description:</b>	<p>Corepoint Scientific™ White Diamond Series Laboratory and Medical Refrigerator  49 Cu. Ft. Capacity  Two swing glass doors, self-closing (energy efficient, double pane vacuum insulated)  Microprocessor temperature controller  Touch Screen Digital Temperature Display with user friendly interface  Alarm Module with battery back-up to support temperature monitoring and alarms  Password-protected control panel  C°/F° switchable  .1° Resolution  Adjustable operating temperature: 1°C to 10°C  Active temperature scale bar display  Resettable 1-point min/max history  Digital 24-hour temperature chart  Temperature data logging  USB port for data transfer and software updates  Real-time operational display clock  4-20mA output  Adjustable alarm volume  Audible and visual high/low temperature alarms  2 probes; 1 in air and 1 in sample bottle  Sensor error alarm &amp; Power failure alarm  Door ajar alarm  Low battery alarm  Alarm validation (audible, visual, and remote alarm)  Alarm mute with ring back  Remote alarm contacts  Ten adjustable shelves (adjustable in ½" increments)  Forced draft air circulation  Keyed door lock  Magnetic door gasket for positive seal  Heavy duty door handles  Swivel casters (3 ½")  LED interior lights are safety shielded and switch controlled  Probe access port (3/8")  White powder-coated steel exterior and interior  Hydrocarbon, natural refrigerant (R290)  Easily cleaned smooth interior  Three year parts and labor warranty, plus an additional two year compressor parts warranty  Exterior dimensions: 54" W x 34 ¾" D x 82 1/2" H with casters; 78 ¼" H without casters  Cycle defrost  115V, 60 Hz, 4.5 Amps, 1/4 HP  ETL listed  ENERGY STAR CERTIFIED  Shipping weight: 538 lbs.</p>