

Product Name:	NSWDUF251JWJ/1
Description:	<p>Corepoint Scientific™ White Diamond Series Ultra Low Temperature Freezer  25 Cu.Ft. Capacity  One solid swing exterior door, right hinged, not reversible  Five insulated inner doors  Microprocessor controller specifications  Touch pad keys and four line display for easy programming  Large LED temperature readout  Audible and visual alarms with silence key for audible alarms  Programmable high and low temperature alarm set points  Programmable access code to limit access to set points and controller configuration  Automatic voltage boost to compensate for low line voltages  Normally open, normally closed alarm contacts and 4-20mA signal available to connect to central alarm system  RS-232 serial output  Battery backup on electronics during power outages  Temperature range: -50°C to -86°C  Four stainless steel shelves, adjustable on ½” centers  Cam action handle for one hand latching of exterior door.  Handle supplied with key lock and door switch for audible and visual alarms  Four heavy duty swivel casters with locks  Vacuum relief port for quick access after door closures  CFC-free refrigerant  Two year parts and labor warranty, plus an additional three year compressor parts warranty  Exterior dimensions: 49"W x 36.3"D x 79.6"H  Interior dimensions: 34.8"W x 23"D x 50"H  2" box capacity: 480 ea.  Upright freezer rack capacity: 30 ea. UF442 racks  230V, 60Hz, 10 amp, HP: ½ x ½  UL/C-UL listed  NEMA plug configuration: NEMA 6-15  Cascade refrigeration system  Shipping weight: 829 lbs.</p>