




COREPOINT
SCIENTIFIC

Exterior
PF051SSSLH/0M


Front View



Side View



Top View



115VAC/60Hz/1PH
NEMA-5-15P

Cu. Ft	Defrost	Door	Int Door	Shelves	W"	D"	H"	Refrigerant	H.P.	Amps	Weight
4.2	Manual	1 Solid LH	0	2 Fixed	23 ³ / ₄	24 ¹ / ₂	33 ³ / ₈	R600a	¹ / ₈	1.5	100 lbs.