

Product Name:	PF051WWWADALH/0M
Description:	<p>White Diamond Series Undercounter Freezer Built-In-LEFT HINGED 4.2 Cu. Ft. Capacity One solid swing door, left hinged Microprocessor temperature controller Digital temperature display Audible and visual high and low temperature alarms Remote alarm contacts 2 probes; 1 in air and 1 in sample bottle Adjustable operating temperature range: -15°C – -25°C Two fixed shelves Keyed door lock Magnetic door gasket for positive seal Recessed door handle Leveling legs Probe access port (3/8") Powder coat exterior High density urethane foam cabinet and door insulation HFC-free refrigerant (R600a), foam insulation and packaging Three year parts and labor warranty, plus an additional two year compressor parts only warranty Pyxis[®], Omnicell[®] and AcuDose RX[®] compatible Exterior dimensions: 23 ¾"W x 24 ½"D x 31 15/16"H Manual defrost 115V, 60 Hz, 1.5 Amps, 1/8 HP ENERGY STAR CERTIFIED ADA compliant UL/C-UL listed Shipping weight: 133 lbs.</p>