

Product Name:	PRF092WWW/0
Description:	<p>Corepoint Scientific™ White Diamond Series Refrigerator &amp; Freezer Combination  9 Cu. Ft. capacity (4.6 Ref./4.2 Freezer)  Two exterior solid swing doors right hinged*  Microprocessor temperature controller, refrigerator and freezer  Digital temperature display with sample simulated probe, refrigerator and freezer  Adjustable operating temperature range: Refrigerator: 1°C – 10°C;  Freezer: -15°C – -25°C  Audible and visual high and low temperature alarms, refrigerator and freezer  Remote alarm contacts, refrigerator and freezer  2 adjustable/1 fixed shelves refrigerator/2 fixed shelves freezer with guard rails on back in Refrigerator to prevent product from being pushed against back wall  Back wall evaporator cover, Refrigerator  Forced draft circulation Refrigerator  Keyed door locks  Magnetic door gasket for positive seal  Recessed door handle, freezer; heavy duty door handle, refrigerator  Leveling legs  Stacking kit  LED interior lights  Two probe access ports (3/8"); one, refrigerator/one, freezer  White powder coated exterior  HFC-free refrigerant (R600a), foam insulation and packaging Three year parts and labor warranty, plus an additional two year compressor parts warranty  Exterior dimensions: 23 ¾"W x 24 ½"D x 66 ¾"H  Manual defrost freezer / cycle defrost refrigerator  115V, 60 Hz, 1.74 Amps Refrigerator, 1.5 Amps Freezer, 1/8 HP  Refrigerator and Freezer  UL/C-UL listed  Shipping weight: 269 lbs.  <b>In stock, with next day shipping</b></p> <p>*Please note this unit requires 2 outlets.  *This product is available left hinged, product number PRF091WWWLH/0</p>