

Product Name:	PRF102WWW/0
Description:	<p>Corepoint Scientific™ White Diamond Series Laboratory Refrigerator &amp; Freezer Combination Freestanding (-20°C)  5.2 cu. ft. capacity (Refrigerator) / 4 cu. ft. capacity (Freezer) One solid swing door (fridge), one solid swing door (freezer)  Microprocessor temperature controller  Digital temperature display  Audible and visual high and low temperature alarms  2 probes; 1 in air and 1 in sample bottle  Remote alarm contacts  Adjustable operating temperature range - Freezer: -15°C – -25°C; Refrigerator: 2°C to 8°C  2 fixed shelves (freezer); 3 adjustable shelves (fridge)  Bottom footed shelf to allow air circulation  Back wall evaporator cover  Keyed door locks  Magnetic door gasket for positive seal  Recessed door handles  Stacking kit  Leveling legs  LED interior light  Probe access port (3/8")  White cabinet and door  High density urethane foam cabinet and door insulation  HFC-free refrigerant (R600a), foam insulation and packaging Two year parts and labor warranty, plus an additional three year compressor parts only warranty  Exterior dimensions: 23 3/4" W x 24"D x 65"H  <i>Note: This unit must have 4" clearance on sides and back for adequate ventilation</i>  Manual defrost (freezer); Cycle defrost (refrigerator)  115V, 60 Hz  1.4 Amps, 1/8 HP (freezer); 1.6 Amps, 1/6 HP (refrigerator)  UL/C-UL listed  Shipping weight: 266 lbs.</p>