

| | |
|---------------|--|
| Product Name: | PRF102WWW/030 |
| Description: | <p>Corepoint Scientific™ White Diamond Series Laboratory Refrigerator & Freezer Combination Freestanding (-30°C) 5.2 cu. ft. capacity (Refrigerator) / 4 cu. ft. capacity (Freezer) One solid swing door (fridge), one solid swing door (freezer) Microprocessor temperature controller Digital temperature display Audible and visual high and low temperature alarms 2 probes; 1 in air and 1 in sample bottle Remote alarm contacts Adjustable operating temperature range- Freezer: -27°C – -33°C; Refrigerator: 2°C to 8°C 3 adjustable shelves (freezer); 3 adjustable shelves (fridge) Bottom footed shelf to allow air circulation Back wall evaporator cover Keyed door locks Magnetic door gasket for positive seal Recessed door handles Stacking kit Leveling legs LED interior light Probe access port (3/8") White cabinet and door High density urethane foam cabinet and door insulation HFC-free refrigerant (R600a) Refrigerator; Natural, Hydrocarbon refrigerant (R290) Freezer Two year parts and labor warranty, plus an additional three year compressor parts only warranty Exterior dimensions: 23 3/4" W x 24"D x 65"H <i>Note: This unit must have 4" clearance on sides and back for adequate ventilation</i> Manual defrost (freezer); Cycle defrost (refrigerator) 115V, 60 Hz 2.2 Amps, 1/4 HP (freezer); 1.3 Amps, 1/6 HP (refrigerator) UL/C-UL listed Shipping weight: 268 lbs.</p> |