

Product Name:	PRF122WWW/0
Description:	<p>Corepoint Scientific™ White Diamond Series Refrigerator & Freezer Combination 12 Cu. Ft. capacity (10.5 Ref./1.7 Freezer)</p> <p>Two exterior solid swing doors (solid door refrigerator/solid door freezer), right hinged, not reversible</p> <p>Microprocessor temperature controller, refrigerator and freezer</p> <p>Digital temperature display with sample simulated probe, refrigerator and freezer</p> <p>Adjustable operating temperature range: Refrigerator: 1°C – 10°C; Freezer: -15°C – -25°C</p> <p>Audible and visual high and low temperature alarms, refrigerator and freezer</p> <p>Remote alarm contacts, refrigerator and freezer</p> <p>Six adjustable shelves with guard rail on back to prevent product from being pushed against back wall, Refrigerator</p> <p>Two adjustable shelves with guard rail on back to prevent product from being pushed against back wall, Freezer</p> <p>Bottom footed shelf to allow air circulation, Refrigerator and Freezer</p> <p>Forced draft circulation Refrigerator</p> <p>Keyed door locks</p> <p>Magnetic door gasket for positive seal</p> <p>Integrated door handles</p> <p>Leveling legs</p> <p>Stacking kit</p> <p>LED interior lights</p> <p>Two probe access ports (3/8"); one, refrigerator/one, freezer</p> <p>White powder coated exterior</p> <p>HFC-free refrigerant (R600a), foam insulation and packaging Three year parts and labor warranty, plus an additional two year compressor parts warranty</p> <p>Exterior dimensions: 23 5/8" W x 26 1/2" D x 80 1/4" H</p> <p>Manual defrost freezer / cycle defrost refrigerator</p> <p>115V, 60 Hz, 2.3 Amps refrigerator, 1 Amp freezer, 1/3 HP refrigerator, 1/8 HP freezer</p> <p>UL/C-UL listed</p> <p>Shipping weight: 289 lbs.</p> <p>In stock, with next day shipping</p> <p>*Please note this unit requires 2 outlets.</p>